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BIBLIOGRAFIA

- 1 Cullen AP. Ozone Depletion and Solar Ultraviolet Radiation: Ocular effects, a United Nations environment programme perspective. *Eye & Contact Lens* 2011;37: 185-190.
- 2 Norval M, Lucas R, Cullen AP, et al. The human health effects of ozone depletion and interactions with climate change. *Photochem Photobiol Sci*, 2011;10:199-225.
- 3 Reiter RJ, Tan DX, Fuentes-Broto L. Melatonin: A multitasking molecule. *Prog Brain Res*, 2010;181:127-151.
- 4 Skene DJ, Arendt J. Human circadian rhythms: Physiological and therapeutic relevance of light and melatonin. *Ann Clin Biochem*, 2006;43:344-353.
- 5 Cancer Council Australia. Slip, Slop, Slap, Seek, and Slide. Available at: www.cancer.org.au/cancersmartlifestyle/SunSmart/Campaignsandevents/SlipSlopSlapSeekSlide.htm. Accessed March 4, 2012.
- 6 Lucas RM, McMichael A, Smith W, et al. Solar Ultraviolet Radiation. Global Burden of Disease from Solar Ultraviolet Radiation. Geneva, Switzerland, World Health Organization, 2006.
- 7 Lucas RM. An epidemiological perspective of ultraviolet exposure — public health concerns. *Eye & Contact Lens*, 2011;37:168-175.
- 8 West SK, Munoz B, Istre J, et al. Mixed lens opacities and subsequent mortality. *Arch Ophthalmol*, 2000;118:393-397.
- 9 Hill D, White V, Marks R, et al. Changes in sunrelated attitudes and behaviours, and reduced sunburn prevalence in a population at high risk of melanoma. *Eur J Cancer Prev*, 1993;2:447-456.
- 10 Tsuchida C, Srivannaboon S, Lim HW. Photoprotection by window glass, automobile glass, and sunglasses. *J Am Acad Dermatol*, 2006;54:845-854.
- 11 Sliney DH. Exposure geometry and spectral environment determine photobiological effects on the human eye. *Photochem Photobiol*, 2005;81:483-489.
- 12 Merriam JC. The concentration of light in the human lens. *Trans Am Ophthalmol Soc*, 1996;94:803-818.
- 13 Javitt JC, Taylor HR. Cataract and latitude. *Doc Ophthalmol*, 1995;88:307-325.
- 14 Diffey BL, Larko O. Clinical climatology. *Photodermatology*, 1984;1:30-37.
- 15 World Health Organization. Global Solar UV Index — A Practical Guide. 2002.
- 16 Sasaki H, Sakamoto Y, Schneider C, Fujita N, Hatsusaka N, Sliney DH, Sasaki K. UV-B Exposure to the Eye Depending on Solar Altitude. *Eye & Contact Lens*, 2011;37: 191-195.
- 17 Glickman RD. Phototoxicity to the retina: Mechanisms of damage. *Int J Toxicol*, 2002;21:473-490.
- 18 Coroneo MT, Muller-Stolzenburg NW, Ho A. Peripheral light focusing by the anterior eye and the ophthalmohelioses. *Ophthalmic Surg*, 1991;22:705-711.
- 19 Podskochy A. Protective role of corneal epithelium against ultraviolet radiation damage. *Acta Ophthalmol Scand*, 2004;82:714-717.
- 20 Abraham AG, Cox C, West S. The differential effect of ultraviolet light exposure on cataract rate across regions of the lens. *Invest Ophthalmol Vis Sci*, 2010;51:3919-3923.
- 21 Coroneo MT, Muller-Stolzenburg NW, Ho A. Peripheral light focusing by the anterior eye and the ophthalmohelioses. *Ophthalmic Surg*, 1991;22:705-711.
- 22 Kwok LS, Daszynski DC, Kuznetsov VA, et al. Peripheral light focusing as a potential mechanism for phakic dysphotopsia and lens phototoxicity. *Ophthalmic Physiol Opt*, 2004;24:119-129.
- 23 Roberts JE. Ultraviolet radiation as a risk factor for cataract and macular degeneration. *Eye & Contact Lens*, 2011;37: 246-249.
- 24 Varma SD, Kovtun S, Hegde KR. Role of ultraviolet irradiation and oxidative stress in cataract formation — medical prevention by nutritional antioxidants and metabolic agonists. *Eye & Contact Lens*, 2011;37:233-245.
- 25 Dillon J, Atherton SJ. Time resolved spectroscopic studies on the intact human lens. *Photochem Photobiol*, 1990;51:465-468.
- 26 Dillon J. Photophysics and photobiology of the eye. *J Photochem Photobiol B Biol*, 1991;10:23-40.
- 27 Glickman RD. Ultraviolet phototoxicity to the retina. *Eye & Contact Lens*, 2011;37:196-205.
- 28 Hu DN, Simon JD, Sarna T. Role of ocular melanin in ophthalmic physiology and pathology. *Photochem Photobiol*, 2008;84:639-644.
- 29 Roberts JE. Ocular phototoxicity. *J Photochem Photobiol B Biol*, 2001;64: 136-143.
- 30 Taylor HR, West S, Munoz B, et al. The longterm effects of visible light on the eye. *Arch Ophthalmol*, 1992;110:99-104.
- 31 Rozanowska M, Jarvis-Evans J, Korytowski W, et al. Blue light-induced reactivity of retinal eye pigment. In vitro generation of oxygen-reactive species. *J Biol Chem*, 1995;270:18825-18830.
- 32 Davies S, Elliott MH, Floor E, et al. Photocytotoxicity of lipofuscin in

- human retinal pigment epithelial cells. *Free Radic Biol Med*, 2001;31:256-265.
- 33 Khachik F, Bernstein PS, Garland DL. Identification of lutein and zeaxanthin oxidation products in human and monkey retinas. *Invest Ophthalmol Vis Sci*, 1997; 38:1802-1811.
- 34 Bernstein PS, Zhao DY, Wintch SW, et al. Resonance Raman measurement of macular carotenoids in normal subjects and in age-related macular degeneration patients. *Ophthalmology*, 2002; 109:1780-1787.
- 35 Chalam KV, Khetpal V, Rusovici R, Balaiya S. A review: role of ultraviolet radiation in age-related macular degeneration. *Eye & Contact Lens*, 2011;37:225-232.
- 36 Cruickshanks KJ, Klein R, Klein BE, et al. Sunlight and the 5-year incidence of early age-related maculopathy: The Beaver Dam eye study. *Arch Ophthalmol*, 2001;119:246-250.
- 37 Taylor HR, Munoz B, West S, et al. Visible light and risk of age-related macular degeneration. *Trans Am Ophthalmol Soc*, 1990;88:163-173.
- 38 Taylor HR, West S, Munoz B, et al. The longterm effects of visible light on the eye. *Arch Ophthalmol*, 1992;110:99-104.
- 39 West SK, Rosenthal FS, Bressler NM, et al. Exposure to sunlight and other risk factors for age related macular degeneration. *Arch Ophthalmol*, 1989;107:875-879.
- 40 Wang JJ, Foran S, Mitchell P. Age-specific prevalence and causes of bilateral and unilateral visual impairment in older Australians: The Blue Mountains Eye study. *Clin Exp Ophthalmol*, 2000;28:268-273.
- 41 Klein R, Klein BE, Knudtson MD, et al. Fifteen-year cumulative incidence of age-related macular degeneration. *Ophthalmology*, 2007;114:253-262.
- 42 Mukesh BN, Dimitrov PN, Leikin S, et al. Five year incidence of age-related maculopathy: Visual impairment project. *Ophthalmology*, 2004;111:1176-1182.
- 43 Wolffsohn J, Eperjesi F, Bartlett H et al. Does Blocking Ultra-Violet Light with Contact Lenses Benefit Eye Health? BCLA Conference, Paper presentation 2012.
- 44 Chandler H. Ultraviolet absorption by contact lenses and the significance on the ocular anterior segment. *Eye & Contact Lens*, 2011;37: 259-266.
- 45 Sliney DH. Intraocular and crystalline lens protection from ultraviolet damage. *Eye & Contact Lens*, 2011;37:250-258.
- 46 Walsh JE, Bergmanson J. Does the eye benefit from wearing ultraviolet blocking contact lenses? *Eye & Contact Lens*, 2011;37:267-272.
- 47 Rosenthal FS, Bakalian AE, Taylor HR. The effect of prescription eyewear on ocular exposure to ultraviolet radiation. *Am J Pub Health*, 1986;76:1216-1220.
- 48 Sasaki K, Sasaki H, Kojima M, et al. Epidemiological studies on UV-related cataract in climatically different countries. *J Epidemiol* 1999;9(Suppl 6): S33-S38.
- 49 Sasaki H, Kawakami Y, Ono M, et al. Localization of cortical cataract in subjects of diverse races and latitude. *Invest Ophthalmol Vis Res*, 2003;44: 4210-4214.
- 50 Hedblom EE. Snowscape eye protection. *Arch Environ Health* 1961;2:685-704.
- 51 Sliney DH. Bright light, ultraviolet radiation Multiple-choice questions — take part at opticianonline.net
- 52 Sliney DH. Eye protective techniques for bright light. *Ophthalmology*, 1983;90:937-944.
- 53 American National Standards Institute (ANSI). American National Standard for Nonprescription Sunglasses and Fashion Eyewear—Requirements. New York, NY, ANSI, Standard Z80.3, 2008.
- 54 British Standards Institution (BSI). Personal Eye Protection — Sunglasses and Sunglare Filters for General Use and Filters for Direct Observation of the Sun. Chiswick, UK, BSI, BS EN-1836, 2005.
- 55 Moore L, Ferreira JT. Ultraviolet (UV) transmittance characteristics of daily disposable and silicone hydrogel contact lenses. *Cont Lens Anterior Eye*, 2006;29:115-122.
- 56 Andley UP, Malome JP, Townsend RR. Inhibition of lens photodamage by UV-absorbing contact lenses. *Invest Ophthalmol Vis Sci*, 2011;52:8330-8341.
- 57 Giblin FJ, Lin L-R, Leverenz VR, Dang L. A class I (Senofilcon A) soft contact lens presents UVB-induced ocular effects, including cataract, in the rabbit in vivo. *Invest Ophthalmol Vis Sci*, 2011;52:3667-3775.

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